

Want Chronic Wounds to Heal Faster?

Use the Power of 100% Pure Collagen.



Collagen Content Matters. The More, The Better.

BIOPAD™ Accelerates Wound Healing

- Type I collagen of equine origin that keeps its native structure¹
- Protects the wound bed from the outer environment, constituting a barrier against exogenous infective agents
- Stimulates the formation of new granulation tissue, the proliferation of fibroblasts, and the deposition of new collagen fibers^{1,2}
- Absorbs wound exudates and can control minor bleeding

BIOPAD™ has the highest collagen content currently on the market (250 mg of collagen)¹ — up to 5 times the standard amount of collagen.¹

Collagen content (%)

BIOPAD	●●●●●●●●●●	100%
ENDOFORM	●●●●●●●●●●	90%
PURACOL	●●●●●●●●●●	88.40%
COLACTIVE	●●●●●●●●●●	83%
PROMOGRAN	●●●●●●●●●●	57.20%

Collagen content (mg)

BIOPAD	████████████████████	~250 mg
COLACTIVE	████████████████████	~200 mg
PROMOGRAN	████████████████████	~60 mg
PURACOL	████████████████████	~50 mg



BIO PAD™

Caring for Wounds

Preparation and Application

■ Preparation:

1. Prepare the wound bed per the facilities wound care protocol and debride the wound when necessary.
2. Hydrate BioPad with sterile saline if the wound has limited exudate. (BioPad will turn into a gel consistency)
3. If the wound has heavy exudate, rinse out the wound bed with saline prior to application.
(Do not hydrate BioPad prior to application)

■ Application:

1. Cut BioPad, if necessary, to fit the size of the wound.
2. Apply directly to the wound bed, covering the entire wound; do not overlap the edges of the wound.
3. Apply a secondary dressing per physician protocol to cover the BioPad.

■ Reapplication:

1. Depending on the amount of exudate, BioPad can be reapplied every day.
2. If the wound has a small amount of exudate, BioPad can be reapplied every 48 hours.
3. BioPad will absorb into the wound bed over time.

■ Indications:

BioPad is indicated for all chronic non-healing wounds that are free of necrotic tissue.

■ Contraindications:

Do not use on patients with a known hypersensitivity to collagen.



PRESENTATION AND SIZES:

BioPad is packaged in a transparent waterproof blister pack and is sterilized by gamma irradiation.

2" x 2" - 3 per box | Order Code: 132622

4" x 4" - 1 per box | Order Code: 132644B

BioPad is covered by Medicare part B - HCPCS code: A6021



AWC Medical Products. All rights reserved.

BIO PAD trademark is owned by or licensed to the AWC Medical Products group of companies. All other trademarks are the property of their respective owners.

©2019

1. Rangaraj A et al (2011): Wounds UK, 7/2: 54-6; 2. Karr JC et al (2011): Adv Skin Wound Care 24: 208-216

BP1-01-042019

For more information contact: www.Awcmed.com 877.250.7391